






OP-LINE IV POLE

Model	Description	LxWxH	Pack	Price €	Picture
		in mm	pz	cad.	
IV Pole in chromed steel					
I-N41121	<p>Material: chromed steel</p> <p>Colour of plastic parts: blue</p> <p>Bottle holder: plastic 4 safety bottle hooks - 2 kg / hook</p> <p>Adjustment: screw safety height adjustment</p> <p>Regulating range: 1.350 - 2.150 mm</p> <p>Tube diameter: 25 / 18 mm</p> <p>Length of outer tube: 1.000 mm</p> <p>Base diameter: 635 mm</p> <p>Base leg profile: 25 x 25 mm</p> <p>Castors: twin castors Ø 50 mm, according to DIN EN 60601 electrically conductive and lockable</p> <p>Total- / base weight (net): appr. 3,5 kg / 1,9 kg</p>	Ø 635 x 1350/2150 h mm		ON DEMAND	
IV Pole in stainless steel					
I-N11121	<p>Material: stainless steel</p> <p>Colour of plastic parts: blue</p> <p>Bottle holder: plastic 4 safety bottle hooks - 2 kg / hook</p> <p>Adjustment: screw safety height adjustment</p> <p>Regulating range: 1.350 - 2.150 mm</p> <p>Tube diameter: 25 / 18 mm</p> <p>Length of outer tube: 1.000 mm</p> <p>Base diameter: 635 mm</p> <p>Base leg profile: 25 x 25 mm</p> <p>Castors: twin castors Ø 50 mm, according to DIN EN 60601 electrically conductive and lockable</p> <p>Total- / base weight (net): appr. 3,5 kg / 1,9 kg</p>	Ø 635 x 1350/2150 h mm		ON DEMAND	
I-N11112	<p>Material: stainless steel</p> <p>Colour of plastic parts: blue</p> <p>Bottle holder: stainless steel to screw in 4 safety bottle hooks - 2 kg / hook</p> <p>Adjustment: single handed safety height adjustment</p> <p>Regulating range: 1.350 - 2.150 mm</p> <p>Tube diameter: 25 / 18 mm</p> <p>Length of outer tube: 1.000 mm</p> <p>Base diameter: 635 mm</p> <p>Base leg profile: 25 x 25 mm</p> <p>Castors: twin castors Ø 50 mm, according to DIN EN 60601 electrically conductive and lockable</p> <p>Total- / base weight (net): appr. 3,5 kg / 1,9 kg</p>	Ø 635 x 1350/2150 h mm		ON DEMAND	

Model	Description	LxWxH	Pack	Price €	Picture
		in mm	pz	cad.	
I-O11112	<p>Effective infusion care by hanging safely IV-bottles and bags</p> <p>Completely stainless - enhanced hygiene and disinfection conditions through smooth, electropolished surfaces</p> <p>Product description</p> <p>Material: stainless steel</p> <p>Colour of plastic parts: blue</p> <p>Bottle holder: stainless steel 4 safety bottle hooks - 2 kg / hook</p> <p>Adjustment: single handed safety height adjustment</p> <p>Regulating range: 1.350 mm - 2.350 mm</p> <p>Tube diameter: 25 / 18 mm</p> <p>Length of outer tube: 1.000 mm</p> <p>Base diameter: 635 mm</p> <p>Base leg profile: 30 x 15 mm</p> <p>Castors: twin castors Ø 50 mm, electrically conductive and lockable</p> <p>Total / base weight (net): appr. 3,6 / 1,9 kg with drip vessel and holder</p>	Ø 635 x 1350/2350 h mm		ON DEMAND	
I-M11121	<p>Material: stainless steel</p> <p>Colour of plastic parts: blue</p> <p>Bottle holder: plastic 4 safety bottle hooks - 2 kg / hook</p> <p>Adjustment: screw safety height adjustment</p> <p>Regulating range: 1.350 mm - 2.150 mm</p> <p>Tube diameter: 25 / 18 mm</p> <p>Length of outer tube: 1.000 mm</p> <p>Base diameter: 635 mm</p> <p>Base leg profile: 30 x 15 mm</p> <p>Castors: twin castors Ø 50 mm, according to DIN EN 60601 electrically conductive and lockable</p> <p>Total- / base weight (net): appr. 3,6 kg / 1,9 kg with drip vessel and holder</p> <p><u>MRI suitability: 1,5 Tesla</u></p>	Ø 635 x 1350/2150 h mm		I-M11121	
I-U02224	<p>Material: stainless steel</p> <p>Bottle holder: stainless steel 4 safety bottle hooks - 5 kg / hook</p> <p>Adjustment: stainless steel screw height adjustment with soft drop system</p> <p>Regulating range: 1.900 - 2.800 mm</p> <p>Tube diameter: 33,7 / 25 mm</p> <p>Length of outer tube: 1.600 mm</p> <p>Base diameter: 800 mm</p> <p>Base leg profile: 45 x 25 mm</p> <p>Castors: twin castors Ø 75 mm, according to DIN EN 60601 electrically conductive and lockable</p> <p>Total / base weight (net): appr. 13,2 kg / 7,5 kg</p>	Ø 800 x 1900/2800 h mm		ON DEMAND	