

## TECHNICAL SHEET



**MAYO TROLLEY**  
**MOD. CSOP**

<b>Concept</b>	Instrument trolley in stainless steel for operating theatre
<b>Structure</b>	In stainless steel AISI 304 Oleodynamic pump in stainless steel Heavy construction, extremely robust and easy to handle
<b>Height</b>	height adjustable by foot pedal from 865 to 1215 h mm
<b>Pedal</b>	The foot pedal can be mounted either on front or back side Pedal with protection rubber
<b>Tray</b>	Removable tray in stainless steel 645x450 mm with rim on four sides
<b>Castors</b>	4 antistatic twin Ø 75 mm wheels with support in polyamide - 2 with brakes Non-marking
<b>Carrying capacity</b>	30 kg
<b>Dimensions</b>	645x450x865/1215 h mm

### FRANCEHOPITAL SAS

Z.I. Ouest - 27 Rue Georges Besse - F-67151 Erstein Cedex / FRANCE - SAS au capital de 40 000€ - RC Strasbourg B 348 139 882

**Italian Branch** 39011 Lana - Bolzano, Zona Industriale 11 - **Logistics Centre** 25012 Malpaga di Calvisano (BS), Via A. Vivaldi 3

[www.francehospital.com](http://www.francehospital.com)