

TECHNICAL SHEET



MAYO TROLLEYS

MOD. CS064 / CS064G

Concept	Instrument trolley in stainless steel for operating theatre
Structure	In stainless steel AISI 304 Oil-pressure piston in special chromed steel Base in tubular frame Ø 25 mm
Height	Adjustable height from 830 to 1170 mm with pedal
Tray	Removable tray in stainless steel 600x400 mm with rim on four sides Mod. CS064 - Fixed tray Mod. CS064G - 360° swivel tray, fixing position by screw
Castors	4 antistatic twin Ø 75 mm wheels - 2 with brakes
Carrying capacity	30 kg
Dimensions	600x400x830/1170 h mm

FRANCEHOPITAL SAS

Z.I. Quest, 27 Rue Georges Besse - F-67151 Erstein Cedex / FRANCE - SAS au capital de 40 000€ - RC Strasbourg B 348 139 882
Italian Branch: I-39011 Lana-BOLZANO, Zona Industriale 11 - **Logistics Centre:** I-25012 Malpaga di Calvisano-BRESCIA, Via A. Vivaldi, 3
www.francehospital.com

TECHNICAL SHEET



MAYO TROLLEYS with Piston in stainless steel

MOD. CS064I / CS064G-I

Concept	Instrument trolley in stainless steel for operating theatre
Structure	In stainless steel AISI 304 Oil-pressure piston in <i>stainless steel</i> Base in tubular frame Ø 25 mm
Height	Adjustable height from 830 to 1170 mm with pedal
Tray	Removable tray in stainless steel 600x400 mm with rim on four sides Mod. CS064I - Fixed tray Mod. CS064G-I - 360° swivel tray, fixing position by screw
Castors	4 antistatic twin Ø 75 mm wheels - 2 with brakes
Carrying capacity	30 kg
Dimensions	600x400x830/1170 h mm