

TECHNICAL SHEET



**MAYO TROLLEYS**

MOD. CS075 / CS075G

<b>Concept</b>	Instrument trolley in stainless steel for operating theatre
<b>Structure</b>	In stainless steel AISI 304 Oil-pressure piston in special chromed steel Base in tubular frame Ø 25 mm
<b>Height</b>	Adjustable height from 830 to 1170 mm with pedal
<b>Tray</b>	Removable tray in stainless steel 700x500 mm with rim on four sides  <b>Mod. CS075</b> - Fixed tray <b>Mod. CS075G</b> - 360° swivel tray, fixing position by screw
<b>Castors</b>	4 antistatic twin Ø 75 mm wheels - 2 with brakes
<b>Carrying capacity</b>	30 kg
<b>Dimensions</b>	700x500x830/1170 h mm

**FRANCEHOPITAL SAS**

Z.I. Quest, 27 Rue Georges Besse - F-67151 Erstein Cedex / FRANCE - SAS au capital de 40 000€ - RC Strasbourg B 348 139 882  
**Italian Branch:** I-39011 Lana-BOLZANO, Zona Industriale 11 - **Logistics Centre:** I-25012 Malpaga di Calvisano-BRESCIA, Via A. Vivaldi, 3  
[www.francehospital.com](http://www.francehospital.com)

TECHNICAL SHEET



**MAYO TROLLEYS**

MOD. CS075I / CS075G-I

<b>Concept</b>	Instrument trolley in stainless steel for operating theatre
<b>Structure</b>	In stainless steel AISI 304 Oil-pressure piston in <i>stainless steel</i> Base in tubular frame Ø 25 mm
<b>Height</b>	Adjustable height from 830 to 1170 mm with pedal
<b>Tray</b>	Removable tray in stainless steel 700x500 mm with rim on four sides  <b>Mod. CS075I</b> - Fixed tray <b>Mod. CS075G-I</b> - 360° swivel tray, fixing position by screw
<b>Castors</b>	4 antistatic twin Ø 75 mm wheels - 2 with brakes
<b>Carrying capacity</b>	30 kg
<b>Dimensions</b>	700x500x830/1170 h mm