



Pi-nero

Multifunction, flexible
and quick-change trolley

The new bi-colour trolley!
A service trolley that can easily transform
into a medication or a therapy trolley (and
many more) thanks to its flexibility and the
available accessories.

Can be equipped with:

- Overbridge and accessory rails
- Accessory tubes
- Snap-on & Snap-off side rails
with the SABU system
- Dozens of accessories: glovebox
holders, dustbins, dispensers, ...

Easy to assemble.

Ships unassembled: knock down
packaging.



Standard features

Posts

extruded aluminium, anodized black

Dimensions

758x532x902 H mm

Shelves

bi-colour ABS black & white

682x450 mm

Rim on 3 sides: 25 H mm

Front rim: 12 H mm

Bottom shelf with built-in bumpers

Median shelf is height-adjustable

Wheels

Single black Ø 100 mm, 2 with brake

>>> all with the original plug-in "CLICK&FIX™" fitting

Weight of trolley

2 shelves: 11.9 kg / 3 shelves: 14.1 kg

Load capacity

Each shelf: 25 kg

Whole trolley: 250 kg

FH-Drawers

610x470 mm drawers – H 76, 156 and 236 mm

(1-2-3 modules)

On plastic coated steel runners

Side accessory rail

UNI 10x25 mm standard, compatible with all accessories, SABU™ fastening system

Overbridge

Original "TWIST™" fastening system

2 Ø 25 tubes are fastened to the posts and can be adjusted in height

10x25 mm extruded aluminium rails are added.

Variable height overbridge – Tubes are shipped inside the posts.

Extraordinary features

Plastic moulding

We use a system that is unique in Italy:

one single mould is used to obtain base, middle and top shelves, all this in **bi-colour**.

The two ABS colours bind seamlessly, and are even stronger than shelves made in a single colour.

Shelves

The retention rim is high enough to contain

objects and prevent them from falling. The 45°

incline allows for perfect cleaning and sanitisation.

Design and production

100% designed and made in Italy.

Customers

Not only hospitals, but also

- Nursing homes
- Clinics and day-care centres
- Dentists
- Veterinary clinics
- Aesthetics clinics
- use in kitchens and gardens

Pi-nero 3 shelves

Service trolley
Cod. PNBW319 & accessories

Dimensions

758x532x902 H mm

Configuration

- 3 shelves
- 1 push handle
- Also available without push handle



Pi-nero 3 shelves

Trolley with 1 FH-Drawer
Cod. PNBW319 & accessories

Dimensions

758x532x902 H mm

Configuration

- 3 shelves
- 1 FH-Drawer
- 1 push handle
- Also available without push handle



Pi-nero 2 shelves

Trolley with 4 FH-Drawers
Cod. PNBW319 & accessories

Dimensions

758x532x902 H mm

Configuration

- 2 shelves
- 2 FH-Drawers of 1 module
- 2 FH-Drawers of 2 modules
- 2 push handles
- Also available without push handles



Pi-nero 3 shelves

Trolley with rails and 1 FH-Drawer
Cod. PNBW319 & accessories

Dimensions

758x532x902 H mm

Configuration

- 3 shelves
- 2 accessory rails
- 1 KOBITE holder for sharp waste containers
- 1 WALLY dustbin
- 1 FH-Drawer of 2 modules
- 1 push handle
- Also available without push handle



Pi-nero 3 shelves + tube

Trolley with accessory tube
Cod. PNBW319 & accessories

Dimensions

758x532x902 H mm

Configuration

- 3 shelves
- 2 push handles
- Accessory tube with glovebox holder with hand sanitiser/soap dispenser
- Also available without push handles



Pi-nero 3 shelves + over

Trolley with overbridge and FH-Drawer
Cod. PNBW319 & accessories

Dimensions

758x532x902 H mm

Configuration

- 3 shelves
- 2 push handles
- 1 FH-Drawer of 2 modules
- Overbridge with 2 accessory rails with 5 tilt-out dispensers with triple horizontal glovebox holder with hand sanitiser/soap dispenser
- Also available without push handles



Pi-nero e-Work

Trolley with overbridge for tablet
Cod. PNBW319 & accessories

Dimensions

758x532x902 H mm

Configuration

- 3 shelves
- 1 push handle
- 1 FH-Drawer of 2 modules
- Overbridge with 2 accessory rails with tablet holder



Pi-nero 2 shelves

Service trolley
Cod. PNBW219

Dimensions

758x532x902 H mm

Configuration

- 2 shelves
- Also available with push handle



medi⁺
HEALTH

Tilaukset: tilaukset@medihealth.fi Tiedustelut: info@medihealth.fi
Medihealth Oy, Haarasuonkuja 5 90240 Oulu, y-tunnus 3131118-6, www.medihealth.fi