

Starsys Medical Carts

Easy to clean advanced polymer material is a marked improvement over traditional metal carts: won't dent, chip, rust, flake or corrode.

Interchangeable 3", 6" and 9" (76, 152 and 230mm) drawers with removable totes can be fully extended for easy access.

Label holders available to quickly identify drawer contents.

Each drawer face has a .5" x 18.375" (12.7 x 466.1mm) polished area allowing for secure adhesion of labels (not supplied by Metro).



The convenient overbridge can be fitted with a wide selection of baskets, shelves, and bins to keep necessary items within reach.

Swing-out Side Storage units increase work surface up to 135%.



Side storage allows easy customization of a wide variety of accessories from trash can and Sharps container to storage bins and shelves.

5" (127mm) polymer caster, 2 directional and 2 total lock casters provide control and stability.



Starsys Carts

Starsys provides a complete system of enclosed carts in a wide variety of heights and system widths. When mobility is a key part of your storage needs, Starsys carts provide a broad selection of unique solutions. This modular system approach allows you to create a cart to meet your exact needs.

Starsys cart widths and depth:

The Starsys cart system is built upon a square module. Multiplying this storage module results in system widths of single 22.75" (578mm), double 42.5" (1080mm) and triple 62.5" (1588mm) widths. The depth on all carts is 24.9" (632mm).

Starsys cart heights:

Starsys carts are available in the following working heights: 28" (711mm), 34" (864mm), 36" (907mm), 39" (983mm), 40" (1016mm), 42" (1060mm), 45" (1136mm), and 48" (1212mm).

Overview:

Starsys carts feature 5" (127mm) swivel casters. The front two casters have a Total-Lock Toe Brake (locks both the wheel and the horn), rear two casters have Swivel-Lock tracking casters.

Starsys carts feature a smooth polymer top/worksurface. If chemical resistance is a concern or heavy equipment is likely to be placed on the top of the cart, Starsys Mobile WorkCenters may be a more appropriate product selection.

An extended mobile base is provided when Heavy-Duty drawers, Active Level shelving or qwikSLOT shelving is specified in a Starsys configuration. The extended mobile base adds an additional 3" (76mm) to the overall depth of the cart. The extended base combines additional counter weights with a deeper footprint to reduce any chance of overbalance when heavy-duty drawers or active level shelving is fully extended.

Indicates Wireless Capabilities

42"

Close-up of system



General Supply with Passive Security
Cat. No. **SXRS40CM1**



Isolation
Cat. No. **SXRSIOS2**



Med Surg
Cat. No. **SXRSMSDRG**



Dressing
Cat. No. **SXRSDRS**

45"



Bedside
Cat. No. **SXRS43CM1**



Bedside
Cat. No. **SXRSBED**



Computer-Ready Bedside
Cat. No. **SXRCOMPBED**



Anesthesia
Cat. No. **SXRSANES** (Key Lock)
Cat. No. **SXRXSANES2** (Electronic Tpad)
Cat. No. **SXRWSANES2** (Wireless Tpad)



Central Line
Cat. No. **SXRXSCLINE** (Electronic Tpad)
Cat. No. **SXRWSCLINE** (Wireless Tpad)



Decentralized Medication
Cat. No. **SXRSDCMED**



Arthroscopic
Cat. No. **SXRSARTH**



I.V. Therapy
Cat. No. **SXRSIV**



Unit Dose Medication
Cat. No. **SXRXSUDOSE** (Electronic Tpad)
Cat. No. **SXRWSUDOSE** (Wireless Tpad)



Trauma
Cat. No. **SXRSTRAMA**

Starsys Carts are built to specifications, therefore, they are non-cancelable, non-returnable.

48"



Biomedical
Cat. No. **SXRSBIO**



Critical Care
Cat. No. **SXRSCCU**



Cast
Cat. No. **SXRSCAST**



Difficult Airway Cart
Cat. No. **SXRSDIFAIR**



SXRS3210L



SXRS3310L

Basic Single Wide, Locking Carts

(in.)	(mm)	Drawer Configuration	Cat. No.
24.88x22.75x41.75	1061x578x1061	1-3", 3-6", 1-9", 0-12"	SXRS1310L
24.88x22.75x41.75	1061x578x1061	3-3", 2-6", 1-9", 0-12"	SXRS3210L
24.88x22.75x41.75	1061x578x1061	5-3", 1-6", 1-9", 0-12"	SXRS5110L
24.88x22.75x44.75	1137x578x1137	1-3", 2-6", 2-9", 0-12"	SXRS1220L
24.88x22.75x44.75	1137x578x1137	4-3", 2-6", 1-9", 0-12"	SXRS4210L
24.88x22.75x44.75	1137x578x1137	6-3", 1-6", 1-9", 0-12"	SXRS6110L
24.88x22.75x47.75	1213x578x1213	0-3", 3-6", 2-9", 0-12"	SXRS0320L
24.88x22.75x47.75	1213x578x1213	2-3", 2-6", 2-9", 0-12"	SXRS2220L
24.88x22.75x47.75	1213x578x1213	3-3", 3-6", 1-9", 0-12"	SXRS3310L

Carts shown in chart have taupe pulls.



SXRS5110L



SXRS1220L



SXRS6110L



SXRS0320L



SXRS2220L

Drawer Pulls (one required per drawer) To specify color for drawer pull, order following catalog numbers for each drawer:



Dark Taupe **SXRDP-TP** White **SXRDP-WH** Cobalt **SXRDP-CB** Red **SXRDP-RE** Orange **SXRDP-OR** Green **SXRDP-GR** Yellow **SXRDP-YL** Slate Blue **SXRDP-BL** Violet **SXRDP-VL** Pink **SXRDP-PK** Black **SXRDP-BK**

Starsys® Carts with Power and Technology provides clinicians with real time access to information with the added benefits of:

- Smaller footprint for tighter spaces and easier maneuvering
- Highly configurable drawers and medication bins
- Built-in Microban® antimicrobial product protection

Cart Options – {46.40}

Cart Bodies with PowerPod™

Description	Cat. No.
40" (1016 mm)H key lock cart body, compact keyboard tray & LiFe PowerPod w/8' (2.4 m) cord	SXRX27K-KLIFE
43" (1092 mm)H key lock cart body, compact keyboard tray & LiFe PowerPod w/8' (2.4 m) cord	SXRX30K-KLIFE
40" (1016 mm)H key lock cart body, compact keyboard tray & LiFe PowerPod, no cord	SXRX27K-KLIFEXP
43" (1092 mm)H key lock cart body, compact keyboard tray & LiFe PowerPod, no cord	SXRX30K-KLIFEXP

*XP models do not include a country specific power cord.
 Note: Powered cart shells include a pullout compact keyboard tray that occupies top 3" (76 mm) of available drawer space.
 Keyboard tray inside dimensions: 13.63"w x16.38"L (346 x417mm).

Power Cord - Country Specific for XP models

Description	Cat. No.
6.5' (2 m) type E/F power cord—Europe, Russia	SXRCORD-EFXP
6.5' (2 m) type G power cord—United Kingdom	SXRCORD-GXP
6.5' (2 m) type I power cord—Australia, China, New Zealand	SXRCORD-IXP
6.5' (2 m) type L power cord— Italy, Chile	SXRCORD-LXP

Drawers

Description	Cat. No.
3" (76mm) Drawer (No Drawer Pull)	SXRS3*
6" (152mm) Drawer (No Drawer Pull)	SXRS6*
9" (229mm) Drawer (No Drawer Pull)	SXRS9*
12" (305mm) Drawer (No Drawer Pull)	SXRS12*

*Starsys standard drawer pull colors available.
 **See previous page for drawer color selection & page 171 for cassette offering

Computer Accessories


Description	Cat. No.
Mouse, Optical, USB	COMP-MSE
Waterproof Medical Antimicrobial Scroll Mouse, Black, USB	COMP-MSEAW
Keyboard, iRocks Ultra X-Slim	COMP-KB15
Cover, Keyboard, iRocks Antimicrobial	COMP-KB15CVR
14.7" L Antimicrobial Cleanable Sealed Keyboard, White, USB	COMP-KB15AC
22" LED DC Monitor 1920 x 1080 with Privacy Screen, Cables	COMP-MON22PS
Corrugated Wire Cord Management 42" (1067mm)	SXRLOOM
Hub, 4 Port, Belkin — Black	SXFLUSBHUB

Computer Mounting Options

Description	Cat. No.
Single Wide Overbridge Computer/Monitor Mounting (Overbridge, Rail, Bracket & Cord Holder)	SXRS505
Double Wide Overbridge Computer/Monitor Mounting (Overbridge, Rail, Bracket & Cord Holder)	SXRD505
Articulating Monitor Arm (2-13 lbs.) (.9-6kg)	SXR310
Articulating Laptop Arm (1-18 lbs.) (.5-8.2kg)	SXR310SEC
Articulating Monitor Arm (2-13 lbs.) (.9-6kg)	SXR318-LD
Articulating Monitor Arm (7.5-25 lbs.) (3.4-11.3kg)	SXR318A-LD



Cart body shown with Drawers, Cassette, Laptop and Arm.



Power Cord Options



PowerPod Interior



Keyboard Tray



SXR310
Articulating Laptop Arm



SXR310SEC + SXR310MSETRAY
Articulating Security Locking Laptop Arm and mouse pad.



SXR318-LD
Articulating Monitor Arm



SXRS505
Overbridge Monitor Rail



Packages include carts and accessories.

Computerized Emergency Department Package
SXRSED-KLIFE



Computerized Bedside Package
SXRSBED-KLIFE



Carts

Cat. No.	Description		
SXR30K-KLIFE	Starsys 43" (1092mm) H, Single Wide Cart, Key Lock, Keyboard Tray and LiFe PowerPod	1	1
SXR420	Advanced Keyless Entry Touchpad	1	1

* -KLIFE carts include 6" (152mm casters) adding 1" to overall cart height.

Drawers

Cat. No.	Description		
SXRS3	3" Drawer-less	2	2
SXRS6	6" Drawer-less	2	2
SXRS9	9" Drawer-less	1	1
SXRDP-BL	Slate Blue Drawer Pulls	5	5

Accessories

Cat. No.	Description		
SXR3LDIV	3" Drawer Divider Long Rail 2pk	2	2
SXR3SDIV	3" Drawer Divider, Short Divider 3pk	4	4
SXR6LDIV	6"/9" Drawer Divider, Long Divider 2pk	1	1
SXR6SDIV	6"/9" Drawer Divider, Short Divider 3pk	2	2
SXR205	Accessory Side Mount Bracket	2	2
FL236	Single Glove Box Holder	1	
SXRBSKT-H	Waste Basket & Holder	1	
SXRPODSLK3	Left Side Swingout Pod with 3 Locking Bins		1
SXRPODSLK1	Left Side Swingout Pod with Top Locking Bin	1	
SXR251	Locking Sharps Cabinet with Container	1	
SXR318	Articulating Monitor Arm (7.5-25 lbs. / 3.4-11.3kg)	1	1

Computer Accessories

Cat. No.	Description		
COMP-MSE	Mouse, Optical, USB	1	1
COMP-KB15	Keyboard, iRocks Ultra X-Slim	1	1
COMP-KB15CVR	Cover, Keyboard, iRocks Antimicrobial	1	1
COMP-AIOTGM24	All-in-one M24 Computer w/ 24" Medical-Grade Touchscreen, i5, 2.9Ghz, 8GB RAM, 240GB SSD	1	1

STARSYS® LOCKS & SECURITY OPTIONS



Lock Type & Style
Information.
Pg. 194-195



Tamper-Evident (Passive Security) Locks

Provides passive security for each drawer individually. Quick visual confirmation of broken seals indicates which drawers have been opened. Includes factory-installed lockbar, drawer locktabs permanently riveted to each drawer, and 100 security seals.

Description	Cat. No.
24"H (610mm) Security Bar, Drawer Tabs & Seals	SXRPSLK-24
27"H (686mm) Security Bar, Drawer Tabs & Seals	SXRPSLK-27
30"H (762mm) Security Bar, Drawer Tabs & Seals	SXRPSLK-30
33"H (839mm) Security Bar, Drawer Tabs & Seals	SXRPSLK-33
36"H (914mm) Security Bar, Drawer Tabs & Seals	SXRPSLK-36



Tamper-Evident Locks

Hinged Push-Button Lockbars

These versatile push button lockbars can be used to lock both drawers and doors. Because the push button mechanism is purely mechanical, no wiring, electronic or batteries are required. Each lockbar includes a backup M300 series keylock, 2 keys, and a passive security locktab (security seals ordered separately: LEC320 = Bag of 100).

For Factory-Assembled Units

Description	Left Mounted Cat. No.	Right Mounted Cat. No.
24"H (610mm) Hinged Lockbar	SXPBL24L	SXPBL24R
27"H (686mm) Hinged Lockbar	SXPBL27L	SXPBL27R
30"H (762mm) Hinged Lockbar	SXPBL30L	SXPBL30R
33"H (839mm) Hinged Lockbar	SXPBL33L	SXPBL33R
36"H (914mm) Hinged Lockbar	SXPBL36L	SXPBL36R

For field retrofit of above key locking bars, add suffix **-KD** to above part numbers.



Hinged Lockbar

PIN Locking

The Starsys PIN Lock is available for Starsys units with polymer drawers or Starsys units with hinged doors. This cost-effective security feature allows users to access the contents without worrying about losing or possessing traditional keys. It provides fast and easy access for up to 19 users and 1 supervisor using unique 4-8 digit codes. PIN lock is low maintenance with 5+ years battery life and hassle free with no software or feature programming making it simple to use and easy to maintain.

Features:

- User PINs are easily programmed to 4-8 digits
- Electronic locking to support up to 20 unique PINs (19 User, 1 Supervisor/Administrator)
- 5+ years battery life, runs on two AAA batteries, easy-to-clean keypad
- Self-lock and manual locking options

Short Door - Single

Cat. No.	Description	Orientation	Hinge	Location
SXRXLHSL	PIN Lock Kit, Self Locking, Left Hinge, Horizontal, for Short Single Door	Horizontal	Left	Left
SXRXLHSR	PIN Lock Kit, Self Locking, Right Hinge, Horizontal, for Short Single Door	Horizontal	Right	Right

Short Door - Double, Center Closing

Cat. No.	Description	Orientation	Hinge	Location
SXRXLSCCBASE	PIN Lock Kit, Self Locking, Base, Double Door, Single Bay	Horizontal	Left/Right	Left/Right

Tall Door - Single *

Cat. No.	Description	Orientation	Hinge	Location
SXRXLVSL	PIN Lock Kit, Self Locking Left Hinge, Vertical, Tall Single Door	Vertical	Left	Left
SXRXLVSR	PIN Lock Kit, Self Locking, Right Hinge, Vertical, Tall Single Door	Vertical	Right	Right

Tall Door - Double Center Closing**

Cat. No.	Description	Orientation	Hinge	Location
SXRXLSCC	PIN Lock Kit, Self Locking, Tall, Double Door, Single Bay	Vertical	Left/Right	Left/Right
SXRCDLOCK-D	PIN Lock Kit, Self Locking, Tall, XD Double Door, Single Bay	Vertical	Left/Right	Left/Right

Overhead Door - Single

Cat. No.	Description	Orientation	Hinge	Location
SXRXLHSL	PIN Lock Kit, Self Locking, Left Hinge, for Single Door	Horizontal	Left	Left
SXRXLHSR	PIN Lock Kit, Self Locking, Right Hinge, for Single Door	Horizontal	Right	Right

Overhead Door - Double, Center Closing

Cat. No.	Description	Orientation	Hinge	Location
SXRXLSCCOH	PIN Lock Kit, Self Locking, Overhead, Double Door, Single Bay	Horizontal	Left/Right	Left/Right



PIN Lock Horizontal & Vertical Mounting

Self-Locking

PIN Locking Continued

Drawers - with Lock Module Above

Cat. No.	Description	Orientation	Hinge Location
SXRCLDWR	PIN Lock Kit, Manual Locking, Bank of Drawers, Doors over Drawers	Horizontal	N/A

Short Door - Single

Cat. No.	Description	Orientation	Hinge Location
SXRCLHML	PIN Lock Kit, Manual Locking, Left Hinge, Single Door	Horizontal	Left
SXRCLHMR	PIN Lock Kit, Manual Locking, Right Hinge, Single Door	Horizontal	Right

Short Door - Double, Center Closing

Cat. No.	Description	Orientation	Hinge Location
SXRCLMCCBASE	PIN Lock Kit, Manual Locking, Base, Double Door, Single Bay	Horizontal	Left/Right

Tall Door - Single*

Cat. No.	Description	Orientation	Hinge Location
SXRCLVML	PIN Lock Kit, Manual Locking Left Hinge, Vertical, Tall Single Door	Vertical	Left
SXRCLVMR	PIN Lock Kit, Manual Locking, Right Hinge, Vertical, Tall Single Door	Vertical	Right

Tall Door - Double, Center Closing**

Cat. No.	Description	Orientation	Hinge Location
SXRCLMCC	PIN Lock Kit, Manual Locking, Tall, Double Door, Single Bay	Vertical	Left/Right
SXRCDLOCK-DML	PIN Lock Kit, Manual Locking, Tall, XD Double Door, Single Bay	Vertical	Left/Right

Overhead Door - Single

Cat. No.	Description	Orientation	Hinge Location
SXRCLOHML	PIN Lock Kit, Manual Locking, Left Hinge, Single Door	Horizontal	Left
SXRCLOHMR	PIN Lock Kit, Manual Locking, Right Hinge, Single Door	Horizontal	Right

Overhead Door - Double, Center Closing

Cat. No.	Description	Orientation	Hinge Location
SXRCLMCCOH	PIN Lock Kit, Manual Locking, Overhead, Double Door, Single Bay	Horizontal	Left/Right

Manual-Locking

* PIN lock for tall single doors apply to standard and XD units.

** PIN lock for tall center closing doors have unique part numbers for standard and for XD units.

Electronic Touchpad Locking Systems

Basic — Our entry level keyless electronic drawer locking system is designed for single-wide carts in heights from 36" to 48" tall.

Features:

- 238 User/Supervisor codes, plus up to 12 Facility Code Access entries
- Auto-relock with software adjustable timeout settings
- Low-Frequency (125 kHz) and High-Frequency (13.56 MHz) Proximity Card Reader Options
- Facility Code Access option allows users to match a generic login based on proximity card data
- Cart audit trail data capability (last 1,500 access events)
- USB connection for cart settings, audit trail reporting, and user updates
- Optional LockView 5Pro Software

Advanced and Advanced Wireless — Full-feature keyless electronic drawer locking system with wireless option is designed for single-wide and double-wide carts in heights from 36" to 48" tall.

Features everything from our Basic system, plus:

- 2988 User/Supervisor codes, plus up to 12 Facility Code Access entries
- Cart audit trail data capability (last 15,000 access events)
- Wireless model allows for cart data exchange over an 802.11 A, B, G, or N (2.4/5.0 GHz) Wi-Fi connection
- Multiple Bay Security

Description	Cat. No.
Basic Keyless Entry Touchpad - Single Bay Locking	SXR412
Advanced Keyless Entry Touchpad - Single Bay Locking	SXR420
Advanced Keyless Entry Touchpad - Double Bay Locking	SXR420-2
Advanced Keyless Entry Touchpad with Wireless - Single Bay Locking	SXR420W
Advanced Keyless Entry Touchpad with Wireless - Double Bay Locking	SXR420W-2
Factory-Installed Low-Frequency (125 kHz) Proximity Card Reader	SXR430PLF*
Factory-Installed High-Frequency (13.56 MHz) iClass/SE/Seos Proximity Card Reader	SXR430PHF*
Standard Low-Frequency (125 kHz) Proximity User Card	SXR-UPROXCRD
Standard Low-Frequency (125 kHz) Proximity Key Ring Tag	SXR-PROXTAG
LockView 5Pro Touchpad Management Software	SXR423

* Proximity card readers are compatible with most industry standard data formats. In some cases, a sample card may need to be submitted to verify compatibility.



SXR412
Basic Keyless Entry Touchpad



SXR420W
Advanced Keyless Entry Touchpad
(Wireless model shown)

Polymer Drawers and Accessories



Removable drawer totes allow for exchange or convenient access outside of cart

Description	Inside Dimensions (Height/Width/Length) in. (mm)	Non-Locking Cat. No.	Locking Cat. No.
Polymer Drawers (Drawer pull required — see selection below)			
3" (76mm) SW Drawer	2.38x16.38x15.88 (60x416x403)	SXRS3	SXRS3-LK
6" (152mm) SW Drawer	5.38x16.25x15.75 (136x412.7x400)	SXRS6	SXRS6-LK
9" (230mm) SW Drawer	8.38x16.25x15.75 (225x412.7x400)	SXRS9	SXRS9-LK
12" (305mm) SW Drawer	11.38x16.25x15.75 (298x412.7x400)	SXRS12	-

Drawer Pulls (one required per drawer) To specify color for drawer pull, order following catalog numbers for each drawer:


Dark Taupe SXRD-TP	White SXRD-WH	Cobalt SXRD-CB	Red SXRD-RE	Orange SXRD-OR	Green SXRD-GR	Yellow SXRD-YL	Slate Blue SXRD-BL	Violet SXRD-VL	Pink SXRD-PK	Black SXRD-BK



SXRLABKIT

Description	Dimensions Height/Length in.	Cat. No.
Drawer Label Holders		
Label Kits (10 Pieces Per Kit)	.75" X 11" (19 x 279mm)	SXRLABKIT

Note: Each drawer face has a .5 x 18.38" polished area allowing for secure adhesion of DYMO and Zebra labels (not supplied by Metro).




SXR3DIV Drawer Divider Kit



SXR6DIV Drawer Divider Kit



SXR3DVR Eggcrate Drawer Divider Kit



Drawer dividers are available to organize your supplies for easy access.

Drawer Accessories

Description	Product Detail (Qty.) in. (mm)	Cat. No.
Drawer Divider Kits		
3" (76mm) Drawer Divider Kit	(2) 16.8" (426.7mm) & (6) 5" (127mm)	SXR3DIV
3" (76mm) Drawer Short Dividers	(3) 5" (127mm)	SXR3SDIV
3" (76mm) Drawer Long Dividers	(2) 16.8" (426.7mm)	SXR3LDIV
Egg Crate Style 3" (76mm) Drawer Divider Kit	(5) 16.8" (426.7mm) & (10) 8" (203.2mm)	SXR3DVR
6" (76mm) Drawer Divider Kit	(2) 16.8" (426.7mm) & (6) 5" (127mm)	SXR6DIV
6" (76mm) Drawer Short Dividers	(3) 5" (127mm)	SXR6SDIV
6" (76mm) Drawer Long Dividers	(2) 16.8" (426.7mm)	SXR6LDIV
Egg Crate Style 6"/9" (152/230mm) Drawer Divider Kit	(5) 16.8" (426.7mm) & (10) 8" (203.2mm)	SXR6DVR

Note: **S** = Short Divider; **L** = Long Divider.

Additional Starsys Totes — Totes provide complete containment for smaller items.



Drawer Totes

Description	Inside Dimensions (Height/Width/Length)		Non-Locking Cat. No.	Locking Cat. No.
	(in.)	(mm)		
3" (76mm) Drawer Tote	2.75x18.25x17.5	70x464x445	SXR3TOTE	SXR3-LKTOTE
6" (152mm) Drawer Tote	5.75x18.25x17.5	146x464x445	SXR6TOTE	SXR6-LKTOTE
9" (229mm) Drawer Tote	8.5x18.25x17.5	216x464x445	SXR9TOTE	SXR9-LKTOTE

Note: For Full Extension Frame and Totes see page 180.

Corner Accessories

Description	Cat. No.
Hospital Grade Power Strip & Cord Wrap	SXR595
I.V. Utility Pole (Attaches directly to left or right rear cart corner)	SXRIV*
Adjustable Sharps Container Bracket Attaches to I.V. Utility Pole	LEC9800
Peel Pouch/Catheter Holder-18" 18" H (457mm) x 6" W (152mm) x 3" D (76mm)	SXR240
Peel Pouch/Catheter Holder-28" 27" H (686mm) x 6" W (152mm) x 3" D (76mm)	SXR243

*Mounting pole included in part number.


SXR243

See page 163 for monitor mounts.

Side Pods — All side pods are 30"H x 19.5"W x 6"L (762 x 483 x 152mm)

Description	Cat. No.
Fixed Side Pod (Non-Locking)	SXRPOD
Fixed Side Pod (with Top-Locking Tilt Bin)*	SXRPODK1
Fixed Side Pod (with 3-Locking Tilt Bins)**	SXRPODK3
Left Swing-out Side Pod (Non-Locking)	SXRPODSL
Left Swing-out Side Pod (with Top-Locking Tilt Bin)*	SXRPODSLK1
Left Swing-out Side Pod (with 3 Locking Tilt Bins)**	SXRPODSLK3
Right Swing-out Side Pod (Non-Locking)	SXRPODSR
Right Swing-out Side Pod (with Top Locking Tilt Bin)*	SXRPODSRK1
Right Swing-out Side Pod (with 3-Locking Tilt Bins)**	SXRPODSRK3

*Top Tray with Dividers and Tilt Bin included.

**Top Tray with Dividers and 3 Tilt Bins included.



Accessorized Pods on Starsys cart

Side Pod Accessories

Description	Dimensions — Inches (millimeters)	Cat. No.
Chart Holder	9.5"H (242mm) x 11.5"W (292mm) x 4.75"L (121mm)	SXRCHRT
Containment Shelf Ledge	1.5"H (38) x 16"L (406mm) x 3.5"D (89mm)	SXRPODLGE
Cup Holder ¹	Cup Sizes: 1.38" (35mm), 2" (51mm), 2.5" (64mm)	SXRCUP
Gas Tank Holder ²	Holds C, D and E Tanks	SXRGAS-H
Glove Box Holder — Single	6.13"H (156mm) x 10.13"W (258mm) x 4"L (102mm)	FL236
Glove Box Holder — Triple	18.25"H (464mm) x 10.13"W (258mm) x 5.5"L (140mm)	FL237
Lockable Sharps Container & Glove Box	18.25"H (464mm) x 13"W (330mm) x 5.5"L (140mm)	SXR251
Lockable Sharps Replacement Containers	20/Pack, 5 Qt. (4.7L) Capacity	SXRMF252
Sharps Container Holder	9.3"H (236mm) x 8"W (203mm) x 3"L (76mm)	SXRSHRPS
Tilt Out Bin	7.5"H (191mm) x 16"L (406mm) x 3.5"D (89mm)	SXRPODBIN
Top Tray (with 3 Dividers) ¹	2.5"H (57mm) x 16.88"W (429mm) x 4.38"L (111mm)	SXRPODTRY
Unit Shelf	4.5"D (115mm) x 16.88"W (429mm) x .63H" (8mm)	SXRPODHLF
Waste Basket & Holder	Holder: 15.5"H (394mm) x 8"W (203mm) x 3"L (76mm)	SXRBSKT-H
Waste Basket Only	12.25"H (311mm) x 11.25"W (286mm) x 8.25"D (210mm)	SXRBSKT

¹Attaches to top of Unit Shelf (SXRPODHLF).

²One or two allowed on Fixed Non-Locking Pods Only.


SXRGAS

SXRPODLGE

SXRPODBIN



SXR205, SXR251



SXR246



SXPULLOUTL



SXPULLOUTR

Side Accessory Bracket

Description	Dimensions — Inches (millimeters)	Cat. No.
Side Accessory Bracket	4.38"H (111mm) x 19.38"W (492mm)	SXR205

Side Accessory Bracket Accessories

Description	Dimensions — Inches (millimeters)	Cat. No.
Chart Holder ¹	9.57"H (242mm) x 11.5"W (292mm) x 4.75"L (121mm)	SXRCHRT
Gas Tank Holder — Single ²	Holds C, D and E Tanks	SXRGAS-H
Glove Box Holder —Single ¹	6.13"H (156mm) x 10.13"W (258mm) x 4"L (102mm)	FL236
Glove Box Holder — Triple ²	18.25"H (464mm) x 10.13"W (258mm) x 4"L (102mm)	FL237
Lockable Sharps Container & Glove Box ¹	18.25"H (464mm) x 13"W (330mm) x 5.5"L (140mm)	SXR251
Lockable Sharps Replacement Containers	20/Pack, 5 Qt. (4.7L) Capacity	SXRMF252
Scope Cabinet (2 Hooks)*	48"H (1220mm) x 9"W (229mm) x 6"L (152mm)	SXR245
Scope Cabinet (4 Hooks)*	48"H (1220mm) x 16"W (406mm) x 6"L (152mm)	SXR246
Sharps Container Holder ¹	9.3"H (236mm) x 8"W (203mm) x 3"L (76mm)	SXRSHRPS
Side Bins (3-Pk) ³	4" (102mm) x 4.5" (114mm) x 18.38" (254mm)	MBP216
Waste Basket & Holder ¹	Holder: 15.5"H (394mm) x 8"W (203mm) x 3"L (76mm)	SXRBSKT-H
Waste Basket Only	12.25"H (311mm) x 11.25"W (286mm) x 8.25"D (210mm)	SXRBSKT

*Accessory Bracket(s) included.
¹Order with (1) SXR205 Accessory Bracket.
²Order with (2) SXR205 Accessory Brackets.
³Order with (3) SXR205 Accessory Brackets


Other Side Accessories

Description	Dimensions — Inches (millimeters)	Cat. No.
Extended Handle	4"H (102mm) x 21.5"W (546mm) x 4.25"L (108mm)	SXPULLOUTL
Extended Handle (Field Retrofit)	4"H (102mm) x 21.5"W (546mm) x 4.25"L (108mm)	SXPULLOUTR
Pullout Writing Surface (Left Side)	17.5"W (445mm) x 15.75"L (400mm); (25 lb. cp.)	SXPULLOUTL
Pullout Writing Surface (Right Side)	17.5"W (445mm) x 15.75"L (400mm); (25 lb. cp.)	SXPULLOUTR

Overbridge and Accessories

Description	Cat. No.
(A) Single-wide Overbridge Assembly with two Hanger Rails**	SXR505
Single-wide Overbridge with Computer/Monitor Mounting Kit	SXR505
Single-wide Overbridge Assembly with no Hanger Rails	SXR505-1000
(B) Hanger Rail (2.19" [55mm] H x 1" [25.4mm] W x 19.5" [495mm] L usable space)	SXR515
Single-wide Overbridge with Top Shelf and Hanger Rail	SXR515
Single-wide Overbridge with 2 Shelves	SXR520
Single-wide Overbridge Shelf (14" [286mm] x 23.25" [590mm])**	SXR560
Double-wide Overbridge Assembly with two Hanger Rails** (for use on double-wide carts only)	SXR570
Double-wide Overbridge with Computer/Monitor Mounting Kit	SXR570
Double-wide Overbridge Assembly with no Hanger Rails	SXR570-1000
Double-wide Hanger Rail (2.19" [55mm] H x 1" [25.4mm] W x 39.5" [1003mm] L usable space)	SXR571
Double-wide Overbridge with Top Shelf and Hanger Rail	SXR571
Double-wide Overbridge with 2 Shelves	SXR571
Double-wide Shelf (43.57" [1106mm] x 10.25" [260mm])**	SXR571
Utility Pole Including Universal Clamp (Attaches to Overbridge)	SXR571
Universal Clamp	SXR571
Short Utility Hook (package of 4) (3.5" [89mm] usable length)	SXR571
Laryngoscope Blade/Specimen Holder	SXR575
(C) Half-Size Metal Utility Bin (5.5" [140mm] H x 5.5" [140mm] W x 5.75" [146mm] L)	SXR581
(D) Full-Size Metal Utility Bin (5.5" [140mm] H x 5.5" [140mm] W x 11.63" [292mm] L)	SXR582
Label/Tape Dispenser (holds up to 2" [50mm] diameter roll with 2.75" [70mm] width of roll space)	SXR583
Wire Supply Basket (7" [178mm] H x 5" [127mm] W x 17" [432mm] L)	SXR585
(E) Utility Bin with Cover (Package of 6)	SXR586
†Sharps Container Bracket	FL584
††Sharps Container Bracket — Holds containers with up to a 41" [1041mm] perimeter	LEC9800
Tilt Bin 3 Units with Brackets (6.75" [172mm] W x 8.13" [206mm] D x 7.5" [191mm] H)	SXR586
Tilt Bin 4 Units with Brackets (5.07" [129mm] W x 7.25" [184mm] D x 6.5" [165mm] H)	SXR586
Tilt Bin 6 Units with Brackets (3.13" [79mm] W x 3.5" [89mm] D x 3.25" [82mm] H)	SXR586

** Maximum weight capacity for overbridge assembly is 40 lbs. (18.1kg).
 *** Maximum weight capacity for overbridge shelf is 40 lbs. (18.1kg).
 † Mounts to universal clamp in addition to hanger rail.
 †† Mounts to I.V. pole or utility pole.




SXR505 Overbridge Assembly (shown with optional accessories)

Starsys® Bin Cassette System

Cassette shown in "Single Bin Access" Mode



4 Level Cassette



Detail of Turn Latch

Cassette shown in "Full Level" Mode



Components — Polymer

Description	No. of Bins Per Level	(in.)	Overall Dimensions Height/Width/Length (mm)	Cat. No.
Cassettes and Bins				
1 Level Cassette Body		5.94x11.75x19.38	151x298x492	SXRCASB1
2 Level Cassette Body		8.94x11.75x19.38	277x298x492	SXRCASB2
3 Level Cassette Body		11.94x11.75x19.38	303x298x492	SXRCASB3
4 Level Cassette Body		14.94x11.75x19.38	379x298x492	SXRCASB4
3" (76mm) wide Cassette Bin	6	2.57x11.94x2.69	65x303x68	SXRBINSB3
4 1/2" (114mm) wide Cassette Bin	4	2.57x11.94x4.07	65x303x103	SXRBINSB4
6" (152mm) wide Cassette Bin	3	2.57x11.94x5.5	65x303x140	SXRBINSB6
8" (203mm) wide Cassette Bin	2	2.57x11.94x8	65x303x203	SXRBINSB8

Note: Bins are shipped with one ID card.

Description	Cat. No.
Bin Accessories	
Cassette Bin Divider Kit (20 dividers)	SXRBINDIV
Cassette Label Kit (includes label strips in Mauve, Slate Blue and Jade for identification of bins)	SXRCASLAB
Cassette Bin Labels (5 sheets of Slate Blue bin labels — 18 of each size per sheet)	SXRBINLAB-BL
Cassette Bin Labels (5 sheets of Jade bin labels — 18 of each size per sheet)	SXRBINLAB-GR
Cassette Bin Labels (5 sheets of Mauve bin labels — 18 of each size per sheet)	SXRBINLAB-RD

SXRBINDIV divides up to 4 bins.

SXRCASLAB includes 3 sheets of labels. One sheet per color. Each sheet contains 18 labels for each size bin. It also includes 18 patient ID bin cards for each size bin and 4 cassette ID labels.



Bins

Medication Cassette Transfer Carts

Description	Overall Dimensions Height/Width/Length (in.) (mm)		Inside Height (in.) (mm)	Approx. Pkd. Wt. (lbs.) (kg)		Cat. No.
Single-Wide, Double-Sided*						
42" (1067mm) Single	41.75x24.88x22.69	1060x632x576	30.25	768	150 68	SXRS40TRAN
45" (1143mm) Single	44.75x24.88x22.69	1136x632x576	33.25	845	150 68	SXRS43TRAN
48" (1219mm) Single	47.75x24.88x22.69	1212x632x576	36.25	921	157 71	SXRS46TRAN
Double-Wide, Double Sided*						
42" (1067mm) Double	41.75x24.88x42.5	1060x632x1080	30.25	768	170 77	SXRD40TRAN
45" (1143mm) Double	44.75x24.88x42.5	1136x632x1080	33.25	845	176 80	SXRD43TRAN
48" (1219mm) Double	47.75x24.88x42.5	1212x632x1080	36.25	921	181 82	SXRD46TRAN

*Cassette storage available on both front and back of cart.



Double-Wide, Double Sided Medication Cassette Transfer Cart
 Cat. No. **SXRD46TRAN**
 (shown with cassettes and bins, sold separately, see above)

Starsys Carts are built to exact customer specifications, therefore, they are non-cancelable, non-returnable.

medi⁺
HEALTH

Tilaukset: tilaukset@medihealth.fi Tiedustelut: info@medihealth.fi
Medihealth Oy, Haarasuonkuja 5 90240 Oulu, y-tunnus 3131118-6, www.medihealth.fi