

Color chart

color chart 1
[atn-07, 08]

			605 snow		1280 porcelain	4635 retro			806 hazelnut
	2795 lemon		614 sunflower	4861 pumpkin		1173 valentine		3105 velvet	4103 fuchsia
4633 apricot		4634 lavanda	2395 blue moon		1277 rain drops		2798 deep sea	2799 universe	
	2389 aquamarine	4035 mermaid		2797 mint	1274 esmerald		4868 pistachio	2394 willow	2794 picnic
2793 chocolate	2390 mocca		3120 silvery		1049 antracit	1661 classy		1208 antistatic	
			10120		20279	07432		20288	10295

color chart 2
[ferrari-stamskin® top]



Stools



SURGEON STOOL

Models: **OP // OPN**



Description

Surgeon stool accomplishes the wide selection of products offered by Novak M for **day surgeries** and ambulatory procedures. Since the stool is comfortable, height adjustable, 360° rotatable, available with or without backrest, it provides a perfectly relaxed working position.

The stool is distinguished by **foot-operated pneumatic height control**, which is an ideal solution, when you need free and clean hands to ensure proper hygiene during surgeries. Instead of the standard **black antistatic padding**, you can also choose from other colors of the top quality water repellent and fire-resistant imitation leather from Novak M color selection.

Five antistatic castors with soft running surface ensure perfect stability during procedures. At the same time the stool is very mobile, so you can effortlessly move it closer to the table and perform your work professionally.

Dimensions

	OP (without backrest)	OPN (with backrest)
base diameter	53 cm	53 cm
seat diameter	36 cm	36 cm
adjustable seat height	48–61 cm	48–61 cm
higher height adjustment	54–73 cm / HOP	54–73 cm / HOP
antistatic castors	Ø 75 mm	Ø 75 mm
safe working load	150 kg	150 kg

Description

Medicine stool is a handy companion that complements your practice. It allows you to work in a comfortable and relaxed position and helps you maintain good posture.

The stool is fully **rotatable** and available with or without a **backrest**. Its height can be easily adjusted to set the appropriate position for your height.

The **star base** allows you to move closer to your patient, while the **round base** provides space for resting your feet. The five **rubber castors** run smoothly and are suitable for all types of floors, allowing you to move from your desk to your patient in no time. **Rubber feet** instead of castors are also available.

For the soft padded seat, you can choose from more than **30 colors of imitation leather** so that it matches the color of your table and other equipment in your practice perfectly.

Equipment and accessories

	round base	star base
base diameter	54 cm	62 cm
seat diameter	36 cm	36 cm
adjustable seat height	46–58 cm / SER-O	51–60 cm / SER-Z
higher height adjustment	53–73 cm / HO	62–82 cm / HZ
castors	Ø 50 mm / SER-KO	Ø 75 mm / SER-KZ
antistatic rubber castors	Ø 50 mm / ACO-K	Ø 75 mm / ACZ-K
rubber feet	SER-N	SER-N
antistatic rubber feet	ACO-N	ACO-N
safe working load	150 kg	150 kg

MEDICINE STOOL

Models:

		without backrest	with backrest
star base	castors	SZ-K	SZN-K
	feet	SZ-N	SZN-N
round base	castors	SO-K	SON-K
	feet	SO-N	SON-N



medi⁺ HEALTH

Medihealth Oy, Haarasuonkuja 5 90240 Oulu, y-tunnus 3131118-6, www.medihealth.fi

Tilaukset: tilaukset@medihealth.fi Tiedustelut: info@medihealth.fi