

# HERA ALUMINUM TREATMENT AND SERVICE CART



The Hera cart is lighter than most steel carts, flexible and easily maneuverable. It is built with an aluminum structure, HPL working surface and dual push handles. The new design prevents the growth of bacteria by using a fully sealed structure, therefore it is very easy to clean

**THE CART WILL MAINTAIN ITS PERFECT AESTHETIC LOOK OVER THE YEARS**

- ✓ **100% Aluminum structure!**
- ✓ **Lightweight**
- ✓ **Flat packed – low carbon footprint**
- ✓ **Sealed structure – no growth of bacteria inside the profiles**
- ✓ **Fully adjustable**

The cart is available in two configurations, a single & double cart as shown below

39" (99 cm) working surface height, 3" casters, and bumpers. 4 & 5" casters are also available (Additional Cost)

Note! The ISO Modules, Dividers and Support Rails are sold separately.



## SINGLE HYGIEIA E CART

Item number: HYG-CRT98-E  
H 105.5 | W 73.9 | D 48.6 cm  
H 41.6" | W 16.4" | D 19.2"  
2 units per box



## DOUBLE HYGIEIA U CART

Item number: HYG-CRT98-UU  
H 105.5 | W 98.5 | D 68.6 cm  
H 41.6" | W 38.8" | D 27"  
1 unit per box



# MODULES, DIVIDERS AVAILABLE



## BELOW ARE AVAILABLE MODULES, DIVIDERS AND ACCESSORIES

STANDARD Storage (thermoplastic polymers):


**CREAM MODULES** ABS




Tray 5cm (2") Deep  
**T054060**



Tray 10cm (4") Deep  
**T104060**



Basket 10cm (4") Deep  
**B104060**



Basket 20cm (8") Deep  
**B204060**


**CLEAR PLASTIC MODULES** Polycarbonate



Tray 5cm (2") Deep  
**T054060PC**



Tray 10cm (4") Deep  
**T104060PC**



Basket 20cm (8") Deep  
**B204060PC**



Tray 20cm (8") Deep  
**T204060PC**

**GREEN MODULES**

INFECTION CONTROL storage (ABS with antibacterial agents):



Tray 10cm (4") Deep  
**T104060ND-AM**



Tray 10cm (4") Deep  
**T104060-AM**




Tray 20cm (8") Deep  
**T204060-AM**

**YELLOW MODULES**

STERILE Storage (can be cleaned and disinfected at high temperature):



Tray 10cm (4") Deep  
**T104060ND HR**



Tray 20cm (8") Deep  
**T204060HR**

## DIVIDERS



long 5cm divider



Short 5cm divider



long 10cm divider



Short 10cm divider



long 20cm divider



Short 20cm divider

CREAM	<b>D0560</b>	<b>D0540</b>	<b>D1060</b>	<b>D1040</b>	<b>D2060</b>	<b>D2040</b>
CLEAR PLASTIC	<b>D0560PC</b>	<b>D0540PC</b>	<b>D1060PC</b>	<b>D1040PC</b>	<b>D2060PC</b>	<b>D2040PC</b>
GREEN	-----	-----	<b>DM1060-AM</b>	<b>DM1040-AM</b>	<b>DM2060-AM</b>	<b>DM2040-AM</b>