

INSTRUMENTAL TABLES



E-01



E-02



E-02




E-03



INSTRUMENTAL TABLES K-4 SERIES MAYO TYPE

E-01KO



E-01KO table is made entirely of stainless steel grade 0H18N9 (AISI 304)

top with deepening, with horizontal 360° rotation, with rotation lock, **hydraulic** height adjustment in the range of 950-1320 mm with foot pedal

base equipped with 75 mm diameter castors with lock

Options for additional charge: removable top for easy disinfection, disassembled by one movement of the lever; double castors of 100 mm diameter with lock

*changing the castors affects the range of the table's height adjustment

Total dimensions:

750 x 500 x 950-1320 mm

[L x S x H]

Top dimensions: 750 x 500 mm

[L x S]



E-01KO with a removable top



Option: a table equipped with double castors of 100 mm diameter



INSTRUMENTAL TABLE E-02 TYPE

Elements of **E-02 tables** made of:

top of stainless steel grade 0H18N9 (AISI 304), with **lifted up edges** or **deepening**, manual height adjustment in the range of 850-1300 mm
base of powder coated steel, equipped with 75 mm diameter castors

Option: table entirely of stainless steel grade 0H18N9 (AISI 304)

Total dimensions:
750 x 500 x 850-1300 mm
[L x S x H]

Top dimensions: 750 x 500 mm
[L x S]



E-02KO with lifted up edges



E-02ST with a deepened top



E-02ST (FR) with a flat top



E-02KO (FR) with a flat top

Total dimensions: 600 x 400 x 900-1500 mm [L x S x H]
Top dimensions: 600 x 400 mm [L x S]

INSTRUMENT TABLE E-03 TYPE

Table E-03 is made entirely of stainless steel grade 0H18N9 (AISI 304)

top with deepening, manual height adjustment by means of a **hand crank** in the range of 850-1200 mm
base equipped with 75 mm diameter castors, all with lock

Total dimensions: 1100 x 600 x 850-1200 mm
[L x S x H]

Top dimensions: 1100 x 600 mm
[L x S]



medi⁺
HEALTH

Tilaukset: tilaukset@medihealth.fi Tiedustelut: info@medihealth.fi
Medihealth Oy, Haarasuonkuja 5 90240 Oulu, y-tunnus 3131118-6, www.medihealth.fi