





*medical tables and chairs*









# BLOODSAMPLE CHAIR

**N**

... for small procedures in laboratory


**POSITIONS**

## DESCRIPTION

- Easy-adjustable arm rests are comfortable and pleasant to touch as well as easy to clean and maintain.
- adjustable arm rests and support handles are attached to chair directly
- Patients settle easily thanks to low chair position and support handles.
- The chair is easily set into trendelenburg position by a shock lever. The lever is pushed by foot so that hands stay free to help the patient.
- The chair can be equipped with legs (STEADY), individually lockable castors (READY) or castors with central locking system (GO)
- thicker seamless padding (6 cm)
- upholstery available in more than 40 colors

## MODEL

BON

## ACCESSORIES

- STU1 – stitched upholstery with memory foam (8 cm)
- STSU1 – stitched upholstery with memory foam (12 cm)
- SIF5 – infusion pole
- PVC-N – PVC cover for leg section
- HCN – head cushion N

